



Liquid Handling Products



cs@irisanalytical.com



+1 706 760 4500



www.irisanalytical.com



The Iris Electronic Pipette Filling Device EF-200 is friendly, efficient and reliable electronic pipette filling instrument. Designed for ease and comfort to suit all types of liquid handling needs in every laboratory. Now, available in a new attractive design.

- Two LED's - Red - low charging - Green - Charging status
- Continuous adjustment of pump speed during operation.
- An universal silicon collet is provided to accept all standards pipettes up to 100 ml. internal knurling on the collet holds the pipette firmly.
- Comfortable grip.
- In-line membrane filter (0.2 micron) prevents **(a)** Aerosol Contamination **(b)** Liquid from entering the equipment by accident.
- Ni-MH 3.6 rechargeable batteries are used in combination with low power consumption circuit and vacuum/pressure pump. allows eight hours of continuous operation.
- Available in attractive colour combinations.
- Battery charger is provided with all models.

Model	Description	Colour	Catalogue No.
EF-200B	Electronic Pipette Filling Device	White body / Grey collet and plunger	4004005-0
EF-200G	Electronic Pipette Filling Device	White body / Violet collet and plunger	4004010-0
EF-200R	Electronic Pipette Filling Device	White body / Red collet and plunger	4004015-0

