



MICROSCOPES



cs@irisanalytical.com



+1 706 760 4500



www.irisanalytical.com



Zoom Stereo Microscope IMS-2000 Series is standard type of microscope where you change magnification using a zoom control knob and have eyepieces and objectives, which can, however, be combined to form a main objective and thus project a separate image of the object in each eye. This allows objects to be viewed three dimensionally.

Suitable for circuit board repair and inspection, SMT work, electronics inspection, coin collecting, gemology and gemstone setting, engraving, repair and inspection of small parts.

Zoom Stereo Microscope IMS-2100 Series

Zoom Stereo Microscope IMS-2100 Series is the standard type of microscope where you change magnification using a zoom control knob. The binocular and trinocular viewing head with wide field eyepieces and has excellent zoom objective to choose 0.7X- 3X and auxiliary objective 2X30 mm, 0.75X120 mm, 0.63X137 mm, 0.5X180 mm objectives and each eye tube has interpupillary adjustment. IMS-2100 Series is capable of providing strong images with a broad field of view, the multipurpose and perfect for basic studies of samples properties. Come with incident and transmitting LED lights, allowing samples to be illuminated for a more comprehensive field of view. Suitable for circuit board repair and inspection, SMT work, electronics inspection, coin collecting, gemology and gemstone setting, engraving, repair and inspection of small parts. Known for its dissecting featuring that allow to examine the huge specimen under 1X -4.5X low zoom magnification objective.



- Zoom objective: 1X-4.5X with optional ZO- 0.7X-3X; Auxiliary objective: 2X30 mm, 0.75X120 mm, 0.63X137 mm, 0.5X180 mm objective
- High point wide field WF10x/23 mm eyepieces
- Maximum magnification of up to 1X-4.5X
- Coarse focusing, focusing range of length 100 mm
- 100V-240V/LED light source with incident and transmitted illumination
- Binocular 45° inclined and a 360° rotatable viewing head with Interpupillary adjusting distance 55-76 mm and trinocular 45° inclined and a 360° rotatable viewing head with Interpupillary adjusting distance 55-76 mm
- Small foot print and light weight
- Easy operating with ergonomic design
- Comes with 1.3MP digital camera and software for quality images
- Optional accessories Photo and Video adaptor

Zoom Stereo Microscope IMS-2200 Series

Zoom Stereo Microscope IMS-2200 Series is the standard type of microscope where you change magnification using a zoom control knob. The binocular viewing head with wide field eyepieces and has excellent zoom objective to choose 0.7X- 3X and auxiliary objective 2X30 mm, 1.5X45 mm, 0.5X165 mm objectives and each eye tube has interpupillary adjustment. IMS-2200 Series is capable of providing strong images with a broad field of view, the multipurpose and perfect for basic studies of samples properties. Come with incident and transmitting LED and Halogen lights, allowing samples to be illuminated for a more comprehensive field of view. Suitable for education, lab research, biology, metallurgy, engineering, chemistry, manufacturing, and in the medical, forensic science and veterinary industries. The microscopes can be used for circuit board repair and inspection, SMT work, electronics inspection, dissection, coin collecting, gemology and gemstone setting, engraving, repair and inspection of small parts. Known for its dissecting featuring that allow to examine the huge specimen under 0.7X-4.5X low zoom magnification objective.



- Zoom objective: 0.7X-4.5X with optional Auxiliary objective: 2X30 mm, 1.5X45 mm, 0.5X165 mm objective
- High point wide field WF10x/20 mm eyepieces
- Maximum magnification of up to 6. 7X-45X
- Coarse focusing, focusing range of length 50 mm
- 12V/15W Halogen and 3W LED light source with incident and transmitted illumination
- Binocular 45° inclined and a 360° rotatable viewing head with Interpupillary adjusting distance 54-76 mm
- Designed with small foot print and light weight
- Easy operating with ergonomic design
- Optional accessories LED ring light, Cold light source, 0.35x, 0.5x C, 1x C-mount

Zoom Stereo Microscope IMS-2300 Series

Zoom Stereo Microscope IMS-2300 Series is the standard type of microscope where you change magnification using a zoom control knob. The binocular and trinocular viewing head with wide field eyepieces and has excellent zoom objective to choose 0.67X-4.5X and auxiliary objective 0.3X and 0.5X objectives and each eye tube has interpupillary adjustment. IMS-2300 Series is capable of providing strong images with a broad field of view, the multipurpose and perfect for basic studies of samples properties. Come with incident and transmitting LED lights, allowing samples to be illuminated for a more comprehensive field of view. Suitable for education, lab research, biology, metallurgy, engineering, chemistry, manufacturing, and in the medical, forensic science and veterinary industries. The microscopes can be used for circuit board repair and inspection, SMT work, electronics inspection, dissection, coin collecting, gemology and gemstone setting, engraving, repair and inspection of small parts. Known for its dissecting featuring that allow to examine the huge specimen under 0.67X-4.5X low zoom magnification objective.



- High point wide field WF10X/22 mm eyepieces
- Maximum magnification of up to 0.67X-4.5X
- Coarse focus knob with tension adjustable, moving range: 106 mm. Focusing arm diameter: 76 mm
- 3W LED light source with incident and transmitted illumination
- Binocular 45° inclined and a 360° rotatable viewing head with Interpupillary adjusting distance 54-76 mm and trinocular 45° inclined and a 360° rotatable viewing head with Interpupillary adjusting distance 55-76 mm
- Small foot print and light weight
- Easy operating with ergonomic design
- Greenough optical system
- accessories adjustable 0.5× C-mount adapter
- Optional accessories Photo Adapter, 1× C-mount adapter

Zoom Stereo Microscope IMS-2400 Series

Zoom Stereo Microscope IMS-2400 Series is the standard type of microscope where you change magnification using a zoom control knob. The binocular and trinocular viewing head with extra wide field eyepieces and has excellent zoom objective to choose 0.8X-5X and auxiliary objective 2X45 mm, 0.5X210 mm objectives and each eye tube has interpupillary adjustment. IMS-2400 Series is capable of providing strong images with a broad field of view, the multipurpose and perfect for basic studies of samples properties. Come with incident and transmitting LED lights, allowing samples to be illuminated for a more comprehensive field of view. Suitable for education, lab research, biology, metallurgy, engineering, chemistry, manufacturing, and in the medical, forensic science and veterinary industries. The microscopes can be used for circuit board repair and inspection, SMT work, electronics inspection, dissection, coin collecting, gemology and gemstone setting, engraving, repair and inspection of small parts. Known for its dissecting featuring that allow to examine the huge specimen under 0.8X-5X low zoom magnification objective.



- Auxiliary 2X45 mm, 0.5X210 mm objective
- Extra wide field EW10X/22 mm eyepieces
- Maximum magnification of up to 0.8X-5X
- Frosted glass plate of length 100 mm, white and black round plate of length 100 mm
- Coarse focus knob with tension adjustable, moving range: 106 mm. Focusing arm diameter: 76 mm
- 100V/240W LED light source with incident and transmitted illumination
- Binocular 45° inclined and a 360° rotatable viewing head with Interpupillary adjusting distance 52-75 mm and trinocular 45° inclined and a 360° rotatable viewing head with Interpupillary adjusting distance 52-75 mm
- Small foot print and light weight
- Easy operating with ergonomic design
- Accessories adjustable 0.5× C-mount adapter
- Optional accessories Photo Adapter

Zoom Stereo Microscope IMS-2500 Series

Zoom Stereo Microscope IMS-2500 Series is the standard type of microscope where you change magnification using a zoom control knob. The binocular and trinocular viewing head with extra wide field eyepieces and has excellent zoom objective to choose 0.8X-5X and auxiliary objective 0.5X220 mm, 0.7X125 mm, 2X45 mm objectives and each eye tube has interpupillary adjustment. IMS-2500 Series is capable of providing strong images with a broad field of view, the multipurpose and perfect for basic studies of samples properties. Come with incident and transmitting LED lights, allowing samples to be illuminated for a more comprehensive field of view. Suitable for education, lab research, biology, metallurgy, engineering, chemistry, manufacturing, and in the medical, forensic science and veterinary industries. The microscopes can be used for circuit board repair and inspection, SMT work, electronics inspection, dissection, coin collecting, gemology and gemstone setting, engraving, repair and inspection of small parts. Known for its dissecting featuring that allow to examine the huge specimen under 0.8X-5X low zoom magnification objective.



- Auxiliary objective 0.5X220 mm, 0.7X125 mm, 2X45 mm objective
- Extra wide field EW10X/22 mm eyepieces
- Maximum magnification of up to 0.8X-5X
- Glass Insert Plate of length 100mm, White and Black round plate of length 100 mm
- Coarse focus knob with tension adjustable, moving range of length 105 mm
- 100V/240W LED light source with incident and transmitted illumination
- Binocular 45° inclined and a 360° rotatable viewing head with Interpupillary adjusting distance 52-75 mm and trinocular 45° inclined and a 360° rotatable viewing head with Interpupillary adjusting distance 52-75 mm
- Small foot print and light weight
- Easy operating with ergonomic design
- Optional accessories photo and video adaptor

