

THERMAL CYCLERS & REAL-TIME PCR SYSTEM



Gradient Thermal Cyclers ITCG-2000 has been used in Clinical detection in hospital, professional, expert. Accommodates 2 different wells including 96×0.2 ml, adaptive to 8-well strip PCR tubes and 96-well plates and this block is completely independently of each other and without interference. Fabricated for easy-to-use, compressed, and accurate thermal cyclers to fit every necessary PCR workflow.

- Operating windows comes with 8 inches TFT vibrant touch graphical display and WinCE
- Compatible with mouse connection with smart touch control as well as access to LAN
- Quicker ramp rates allow more work to be done on a workday
- Simple to program unit is compact in design and built to perform
- Bluetooth printer for easy operation
- Meet the need for another testing hot lid temperature and hot lid work mode available
- 4 Temperature Control modules
- Great reproducibility of PCR reaction
- Achieve traceability of experimental result
- Portable and User-friendly
- Faster speeds
- No condensation after overnight cooling at $\geq 4^{\circ}\text{C}$

