

WEIGHING SCALES/BALANCES



Top-loading Balances IBL-3000S Series feature weighing capacities of 100 g, 200 g, 300 g, 500 g, 600 g, 1000 g, 2000 g and readabilities of 0.01g. Some applications include weighing, counting, check weighing, weighing under unstable conditions, peek hold and percentage (%).

Iris top-loading balances are designed for routine weighing applications in your workplace, offering accuracy and repeatability during essential weighing applications in education, laboratory and, industrial sector. It delivers competitive performance at an economical price, featuring RS232 connectivity for easy communication, and a backlit display and a simple interface for uncomplicated operation, these balances are perfectly designed for your workplace.

Standard features include Basic functionality for routine weighing applications. Equipped with three essential weighing modes, RS232 connectivity for data transfer and storage and internal calibration (in selected models) making them ideal for routine weighing applications.

Designed for uncomplicated operation with easy-to-use display and interface Equipped with an easy-to-read, bright backlit display and a simple user interface, they are incredibly easy to operate, with almost no training required.

Smart design and durable construction Iris top pan balances have a small footprint which saves on benchtop space while providing a large weighing surface. They are durably constructed, and feature a stainless-steel pan to withstand day-to-day use in the workplace.

- Standard RS-232 connectivity for easy data transfer
- Essential weighing modes
- Easy-to-read, bright backlit display and a simple user interface
- LCD display for showing the active functions
- High readability 0.01g
- Fast stabilization and response time
- 130 mm and round weighing pan
- Dry Battery Box, Weight(1g~500g balance)
- Wind Shield, under hook, rechargeable battery, printer(dot-matrix printer, thermal printer)(optional)

Model	IBL-3050S		IBL-3150S		IBL-3250S		IBL-3350S		IBL-3450S		IBL-3550S		IBL-3650S	
Capacity	100 g		200 g		300 g		500 g		600 g		1000 g		2000 g	
Display	LCD													
Readability	0.01g													
Pan Dimension	130 mm													
Measuring Parameters	g, ct, oz, lb, GN, %, pcs													
Interface	RS232													
Certification														
Packing Dimension	330x280x150 mm													
Gross Weight	3 kg													
Electrical Requirement	220V	110V	220V	110V	220V	110V	220V	110V	220V	110V	220V	110V	220V	110V
Catalogue No.	6801200502	6801200501	6801201002	6801201001	6801201502	6801201501	6801202002	6801202001	6801202502	6801202501	6801203002	6801203001	6801203502	6801203501

