

WEIGHING SCALES/BALANCES



Moisture Balance IBMH-4000 Series measures moisture content with increased precision and efficiency with halogen-infrared heat source features, PT100 sensor to ensure process repeatability and wide data storage capacity along with touch panel to display time, temperature and weighing results all at a time. Application area including biology, food, agriculture, pharmacy, textile, dye and research industries to resolve moisture content in a product. A single 400 W halogen lamp heats samples at 1 °C selectable rises. Equipped with color keys to enable quick recognition. Temperature range 180°C and have capacity to weigh 50 g with readability of 1 mg, 5 mg, 10 mg.

Iris moisture balance is perfect for measuring moisture in most products including grains, pharmaceuticals, soils, sludge and chemicals. The amount of moisture in a product can greatly affect perceived quality. This measurement is key in quality control process of material as it qualifies quality or purity of a product.

There are different types of heating elements to be considered: ceramic, halogen and quartz elements. Each element offers heating properties so you should consider your application before choosing the appropriate one. Our tech support team can help you to select the relevant model.

The principle of moisture balance is to measure the original weight of a product before the drying process. It automatically and simultaneously dries and weighs solid samples for the determination of moisture content. It provides a continuous direct readout for both weight and percentage moisture loss through the entire cycle

In research work, manufacturing and quality control applications it is essential to be able to determine moisture content simply and reliably. This moisture balances do precisely that, providing moisture content data as a percentage of the initial mass or the residual dry mass.

Model	IBMH-4300	IBMH-4600	IBMH-4900
Capacity	50 g	50 g	50 g
Moisture Range			
Temperature Range	180?		
Dry Residual			
Resolution	1 mg	5 mg	10 mg
Display	Digital		
Heat Source	Halogen lamp		
Certification			
Time Range	3- 4 second		
Power Consumption			
Dimension	300x160x140 mm		
Weight			
Electrical Requirement of Balance			
Catalogue No.	68015505-04	68015510-04	68015515-04

