



STEAM STERILIZERS/AUTOClaves



cs@irisanalytical.com



+1 706 760 4500



www.irisanalytical.com

Tabletop Autoclave ITA-200 Series

The ITA-200 front loading 'Benchtop' autoclave range is available in four sizes. These units require no plumbing, and only a standard single phase power supply. A 5.7-inch colour touchscreen with integrated USB data port. Electro polished stainless steel chamber. Simple pre-post vacuum option available with adjustable time and temperature range. The ITA-200 is available in four sizes, 12l, 16l, 18l, and 23liters and is one of the most energy efficient autoclaves available. The pressure vessel design ensures the absolute minimum amount of water is used during a sterilization cycle, which not only means comparably low power consumption, but also very quick heating and cooling. A benchtop autoclave has a variety of uses, being suitable for sterilizing laboratory glassware, plastic receptacles and containers, petri dishes, un-wrapped instruments, media, and occasional discard loads.

In addition to laboratory use, it is widely used in scientific research, process manufacturing, veterinary and medical healthcare industries.

- LCD display panel
- Adopt Motorola sensor for intelligent and precise control as well as fault auto-detection technique
- Precise temperature control ($\pm 0.5^{\circ}\text{C}$) during sterilization process
- Built-in program to clean the pipelines automatically and avoid the blockage
- The LCD can display the sterilization data and relevant curve clearly
- Built-in water quality, water shortage and water overrun monitoring features
- Safety value, double lock, mechanical lock, air lock and multi-level protection to ensure the pressure is safe
- Optional for printer, built-in USB module to output the sterilization data via USB Flash.
- Widely applied in orthopaedics, ophthalmology and beauty etc. industries.
- Micro-computer control, imported Motorola sensor, high accusation and precision.
- During the sterilizing procedure, the range of temperature is ± 0.5 .
- Automatic cleaning program in the steamer, avoid the occlusion.
- Display the sterilizing procedure date and curve.
- Water quality monitoring, indicate the water tank condition, automatic notice when used water tank is full.
- Safety valve, double lock, mechanical lock, air lock, multi pressure safe protection device.
- Printer is optional, record the sterilizing procedure with USB.
- Applied with orthopaedics, ophthalmology, beauty industries and other area.
- Built-in USB interface, can store up to 2000 data cycles
- Sterilisation data can be transferred to computer/printer
- New Tabletop Water tank design
- Pressure, mechanical, electronic door lock system, completely cut off hidden danger.
- Double temperature monitoring system, analyze and compare the different temperature of the tank, to guarantee and control every heating stage accurately.
- Adopt high analytical pressure sensor, the pressure can accurate to 0.001 bar, in the mean time we add sampling system to improve the accuracy of information.
- Circular heating reaches minimum temperature difference. Pressure-prevention and corruption-resistance makes it highly reliable. Alloy heating loop protects against damage by no-water heating.
- High efficiency double pump, low noise, vacuum ability can reach -0.92 bar, double head ensure the balance during working and stable vacuum state in the chamber.
- Raise 10 times temperature controlling precision system, can accurate to 0.1°C .
- Stable ULKA water pump
- Pass by all kinds of bacteria test, come with B&D testing system, vacuum testing system, Helix testing system, can examine every parameter index, to guarantee the validity of sterilization.
- Original fuzzy control system driven by steam generator adopts OWM(pulse with modulation), able to discriminate the fine differences in temperature allowances and to adjust the modulated waves correspondingly thus making temperature rise at the same speed each time
- High quality CEME imported valve

Model	ITA-212		ITA-216		ITA -218		ITA-223	
Chamber Volume	12 L		16 L		18 L		23 L	
Insulation Class	class B							
Control Accuracy	Temperature 0.1 ℃, Pressure 0.001 bar							
Sterilization Level	Class B							
Sterilization Temperature	121?,134?							
Special Sterilization	HIV/AIDS, hepatitis b (HBV), mad cow virus and spore.							
Drying Procedure	A strong negative vacuum drying							
Operation Interface	Adopt the operation of the new interface LCD screen can display temperature and pressure time running state, fault alarm information, etc							
Test system	Bowie & Dick Steam penetration test, Vacuum test, Helix test							
Control Accuracy	Temperature 0.1℃, pressure 0.1 K pa.							
Sterilization Data	External USB connection , (printer is optional)							
Multiple Safety	Interlock protection device: Safety valve, Pressure interlock, door alarm, automatic breakdown check. Over temperature and pressure protected.							
Water System	Internal open water tank, easy to clear.							
Sterilization Configuration	Activities tray rack, match 3 pallets							
Display	LCD display							
Working Temperature	134 ℃							
Working Pressure	2.1 bar							
Warm Up Time	0 to 7 minutes							
Pump down Time	10 minutes							
Air Intake Time	4 minutes							
Working Time	6 minutes							
Drying Time	18 minutes							
Quick Drying Time	9 minutes							
Total Time	40 minutes							
Security	Safety valve, pressure system, mechanical and electric lock. Self-diagnosis, pressure-protection device, heating-exceeding							
Test Programmers	Helix/BD Test; Vacuum Test							
Drying Progress	3-times pulsation vacuum, temperature less than 0.2% after finish							
Interface	USB/External Printer							
Power Consumption	1600 W		1700 W		1800 W		2000 W	
External Dimensions	625x425x456 mm		625x425x456 mm		625x425 x456 mm		625x425 x456 mm	
Chamber Dimensions	200x360 mm		230x360 mm		250x360 mm		250x450 mm	
External Dimensions	616x460x420 mm		616x460x420 mm		616x460x420 mm		616x460x420 mm	
Net Weight	-		-		-		-	
Shipping Dimensions	655x545x450 mm		655x545x450 mm		655x545x450 mm		750x545x450 mm	
Gross Weight	-		-		-		-	
Electrical Requirements	110V/60Hz	220V/50Hz	110V/60Hz	220V/50Hz	110V/60Hz	220V/50Hz	110V/60Hz	220V/50Hz
Catalogue No.	580050501	580050502	580051001	580051002	580051501	580051502	580052001	580052002

