



STEAM STERILIZERS/AUTOClaves



cs@irisanalytical.com



+1 706 760 4500



www.irisanalytical.com

Benchtap Autoclave/Steam Sterilisers ITA-700 Series provide the broadest range of steam sterilizers and has a solution available to maximize department productivity and optimize workflow. It can be used to Inactivated AIDS (HV), hepatitis B (HBV) mad cow virus and spores. The pressure door is auto locking device provide easy option for such a cost effective autoclave. Design with water supply and detection system that is used to clean disinfection with functions of water quality detection and debris filter. Mainly used in Dental clinic, medical, laboratory, Beauty Shop applications. Equipped with excellent features contain relief valve, pneumatic locking, and breakdown self-checking system with LCD display of pressure and temperature.

- Versatile and compact design structure
- Corrosion resistant
- Temperature stability
- European Standards (EN 13060)
- Build in LCD display
- Automatic door locks
- Smart and advance operating interface
- Available with pipe connection
- Pipe materials is metal constructed provide excellent durability and consistency, bear at high temperature
- Automatic door lock system
- Water supply and detection system

Model	ITA-745
Chamber	45L
Sterilization Temperature	121°C, 134°C
Tray	2-layer movable tray
Level	European Class B Standard
Safety	Safety valve, pneumatic lock, fault self-checking system
Test	Bowie & Dick steam penetration test, helix spiral test, Vacuum test, water quality testing
Drying program	Powerful vacuum drying, the remaining humidity of the instrument is <0.2%
Display	LCD
Applications	Inactivated AIDS (HV), hepatitis B (HBV) mad cow virus and spores
Interface	printer and USB
Power Consumption	3500 VA
External Dimension	610×840×580 mm
Internal Dimension	319×620 mm
Net Weight	100 kg
Electrical Requirements	AC 220V/50Hz
Catalogue Number	580300502

