



GLASSWASHERS and DISINFECTORS



cs@irisanalytical.com



+1 706 760 4500



www.irisanalytical.com



Laboratory Glassware Washers IGW-1000 Series suitable for washing and drying narrow-neck with HEPA filter, general purpose glassware washer and has specialized features for contamination-sensitive research especially for laboratory and medical applications. The IGW-1000 is decontaminating equipment that worked on moist heat to remove contamination and bacteria with a large wash cabinet and short programme cycles. The key features include stainless steel interior, exterior, racks and inserts withstand heavy use and provide protection against corrosion resistance materials and offer internal components peristaltic pumps and seals laboratory-grade for excellent durability of instruments.

- With default and standard programs with user defined and custom programming
- Washing chamber and internal door are made up to 316L and 304L stainless steel
- Provide over temperature and leakage protection
- Large capacity due to optimised chamber dimensions
- Versatility and reliable
- Use peristaltic pump for liquid detergent and neutralizer
- Independent air heater, HEPA filter and drying process
- Quick and effective drying system
- Offer fault diagnosis, electronic lock, automatic power off function
- Comes with rotatory knob and 4 castors
- Suitable for cleaning laboratory glass vessels, ceramic, metal and plastic materials

Accessories

- RS 232
- SD card
- Printer
- Pure booster pump

Model	IGW-1200		IGW-1400	
Function	Washing		Washing + Drying with HEPA filter	
Capacity	176 L		176 L	
Circulating System	0 to 400 L/min		0 to 400 L/min	
Materials	Internal Panel: 316L, External Panel: 304L		Internal Panel: 316L, External Panel: 304L	
Peristaltic pump	2 to 4 pcs		2 to 4 pcs	
Programs	8 Standard, 2 user defined		8 Standard, 2 user defined	
Function	Fault diagnose, electronic lock, automatic power off		Fault diagnose, electronic lock, automatic power off	
Interface	RS232, SD card, printer, pure booster pump		RS232, SD card, printer, pure booster pump	
Heating Power	4 KW/9 KW		4 KW/9 KW	
Dimension	950×605×750 mm		950×925×750 mm	
Electrical Requirement	220V/50Hz	380V/50Hz	220V/50Hz	380V/50Hz
Catalogue No.	370050502	370050503	370051002	370051002

