



GLASSWASHERS and DISINFECTORS



cs@irisanalytical.com



+1 706 760 4500



www.irisanalytical.com



Laboratory Glassware Washers IGW-4000 closed and compact sophisticated model and easy to use. The chamber capacity is up to 308 L and perfect solution for both small and large items as it has countertop installation as well as under counter installation. The IGW-4000 is decontaminating equipment that worked on moist heat to remove contamination and bacteria with a large wash cabinet and short programme cycles. The key features include stainless steel interior, exterior, racks and inserts withstand heavy use and provide protection against corrosion resistance materials and offer internal components peristaltic pumps and seals laboratory-grade for excellent durability of instruments.

- With default and standard programs with user defined and custom programming
- Washing chamber and internal door are made up to 316L and 304L stainless steel
- Provide over temperature and leakage protection
- Automatic safety lock
- Large capacity due to optimised chamber dimensions
- Versatility and reliable
- Use peristaltic pump for liquid detergent and neutralizer
- Quick and effective drying system
- Suitable for cleaning laboratory sampling tubes, measuring tubes, pipette in bulk quantity
- Comes with automatic door function
- Integrated with PT1000 sensor
- Offer interchangeable cleaning frame and baskets
- Comes with automatic sliding glass door
- Provide steam generator
- Comes with countertop and under counter installation
- Provide conductivity monitoring system

Accessories

- RS 232
- Pure booster pump

Model	IGW-4000	
Capacity	308 L	
Materials	? 95°C	
Peristaltic pump	Internal Panel: 316L, External Panel: 304L	
Cleaning Quantity	PT1000	
Programs	Sampling Tube: 672 pcs, Measuring Tube 100 (ml): 126 pcs, Pipette: 238 pcs	
Function	36 Default, 100+Custom	
Heating Power	RS232, pure booster pump	
Total Power	3/8/18 KW	
Dimension	5/11/20 KW	
Weight	1350×925×760 mm	
Electrical Requirement	220V/50Hz	380V/50Hz
Catalogue No.	370200502	370200503

