



GLASSWASHERS and DISINFECTORS



cs@irisanalytical.com



+1 706 760 4500



www.irisanalytical.com



Laboratory Glassware Washers IGW-5000 is a double glass sliding door large glassware washer for cleaning both small and large items and increase workflow and output as user can load 5 independent racks especially for laboratory and medical applications. The IGW-5000 is decontaminating equipment that worked on moist heat to remove contamination and bacteria with a large wash cabinet and short programme cycles. The key features include stainless steel interior, exterior, racks and inserts withstand heavy use and provide protection against corrosion resistance materials and offer internal components peristaltic pumps and seals laboratory-grade for excellent durability of instruments.

- With default and standard programs with user defined and custom programming
- Washing chamber and internal door are made up to 316L and 304L stainless steel
- Automatic safety lock
- Large capacity due to optimised chamber dimensions
- Versatility and reliable
- Use peristaltic pump for liquid detergent and neutralizer
- Quick and effective drying system
- Comes with automatic door function
- Integrated with PT1000 sensor
- Independent air heater, HEPA filter and drying process
- Comes with automatic sliding double glass door
- Use 5 level loading baskets
- Provide dosing detection
- With spray arm speed monitoring and basket self-recognise function

Model	IGW-5000
Capacity	7480 L
Circulating System	1300L/min.
Basket Loading	5 Levels
Materials	Internal Panel: 316L, External Panel: 304L
Programs	36 Default, 100 user defined
Heating Power	9/18/27 KW
Total Power	12/21/30 KW
Environment Temperature	5 to 40°C
Dimension	2340×1260×1050 mm
Electrical Requirement	380V/50Hz
Catalogue No.	370250502

