

POLARIMETERS & REFRACTOMETERS





Semi-Automatic Refractometer used in the manufacturing center and production line seeking to maintain a high level of quality and adherence to stringent quality procedures, allow accurate reading of data and still require some manual adjustment to handle this equipment and measuring data rely on user skill.

Semi-Automatic Refractometer ISAR-1000 Series

Semi-Automatic Refractometer ISAR-1000 Series used in the manufacturing center and production line seeking to maintain a high level of quality and adherence to stringent quality procedures, allow accurate reading of data and still require some manual adjustment to handle this equipment and measuring data rely on user skill. Provide accurate readout of the refractive index, %brix value, and temperature that display on LCD touch display. The ISAR-1000 series provide refractive index range, up to 1.30000 to 1.70000, resolution 0.00001 and brix range, up to 100 %, resolution between 0.1 and 0.01%, with temperature range 0 to 50°C and -20 to -120°C. These refractometer normally have RS232 and USB interface with a PC and are the instruments of choice for high-throughput labs seeking both accuracy and data integrity.



- Used in petroleum, chemical, pharmaceutical, sugar refining and food industries, as well as in related colleges, universities and scientific research institutions
- Suitable to measure Refractive index and Brix value
- Offer measuring range of 1.30000 to 1.70000 nD and 0 to 100% Brix
- Easy and portable handling
- Build-in digital display to control and show all parameters
- Provide easy calibration
- High quality strong housing materials

