

CENTRIFUGES



cs@irisanalytical.com



+1 706 760 4500



www.irisanalytical.com

Portable Centrifuge IPC-100 Series are ideal for the separation of serum plasma, urea, blood samples and other routine applications in hospital and research laboratories. The IPC-100 are most powerful and competitive centrifuge in its class. it can be placed with 6 of 15ml conical tubes or glass tubes, 8 of 15ml conical tubes or glass tubes conical tubes which with screw caps The speed can be 100- 5000RPM. the LCD display with ability to change be between RPM and RCF.

Model	IPC-110		IPC-120		IPC-130	
Max. capacity	15mlX6		15mlX8		5mlX12	
Max. speed rpm (Angle rotor)	5000 rpm		5000 rpm		5000 rpm	
Speed increment	100 rpm		100 rpm		100 rpm	
Max. RCF[xg]	3074xg		3074xg		3074xg	
Noise level at max. speed(approximately)	?60dB(A)		?60dB(A)		?60dB(A)	
Max. Acceleration time	20-80 seconds		20-80 seconds		20-80 seconds	
Max. Deceleration time	20-80 seconds		20-80 seconds		20-80 seconds	
Power consumption drive	40W		40W		40W	
Electrical Requirements	110V/60Hz	220V/50Hz	220V/50Hz	220V/50Hz	110V/60Hz	220V/50Hz
Catalogue No.	140200501	140200502	140201001	140201002	140201501	140201502

