

CENTRIFUGES



cs@irisanalytical.com



+1 706 760 4500



www.irisanalytical.com

Portable Centrifuge IPC-200 Series ideal for the separation of serum plasma. urea, blood samples and other routine applications in hospital and research laboratories. A small footprint and quiet operation makes it perfect for desktop use or quick transport between benchtops. The IPC-200 is the most powerful and competitive centrifuge in its class. It can be placed with 8 of 5ml or 4 of 100ml conical tubes or glass tubes conical tubes which with screw caps The speed can be set from 1000 to 14500rpm and extremely compact with exceptionally small footprint, spin up to 14,500 rpm, (14,100 x g).

Model	IPC-210		IPC-210 Plus	
Max. capacity	400ml(5ml×8/100ml×4)		50ml(5ml×10)	
Max. speed rpm (Angle rotor)	13,400 rpm		14,500 rpm	
Speed increment	1000 rpm		1000 rpm	
Max. RCF[×g]	12,100 x g		14,100 x g	
Noise level at max. speed(approximately)	?60dB(A)		?60dB(A)	
Max. Acceleration time	<13s		<13s	
Max. Deceleration time	<13s		<13s	
Power consumption drive	100W		105W	
Weight including rotor (Kg)	4.3kg		4.3kg	
Dimensions (D×W×H) mm	260 x 260 x 150 mm		260 x 260 x 150 mm	
Electrical Requirements	110V/50Hz	220V/50Hz	110V/60Hz	220V/50Hz
Catalogue No.	140250501	140250502	140251001	140251002

