

CENTRIFUGES



cs@irisanalytical.com



+1 706 760 4500



www.irisanalytical.com



Blood type ID Card Centrifuge builds with multi-color cards rotors accessories. Suitable for blood group testing, micro column gel immune detection. Driven by brushless motor. Auto-balancing centrifuge for blood bank is particularly designed for the quick centrifuge in blood bank.

Blood type ID Card Centrifuge IIDC-1000

Blood type ID Card Centrifuge IIDC-1000 builds with multi-color cards rotors accessories. Suitable for blood group testing, microcolumn gel immune detection. Spins sample volumes up to 24 cards micro gel rotor with a maximum speed of 4000 rpm. Driven by a brushless motor. Auto-balancing centrifuge for blood bank is particularly designed for the quick centrifuge in blood bank with storages up to 9 program records for users. For safety, we introduced an electronic door lock function and, imbalance detection, over-speed protection, inching function, short-spin function.



- Specifically for blood group test, microcolumn gel immune detection
- The machine adopts brushless motor driving, micro-computer control
- Electric lid lock, super speed, and self-test protection
- Low noise, without carbon dust pollution
- Professional program design Can be run directly Without tedious parameter settings
- Maximum capacity of up to 24 cards micro gel rotor

Blood type ID Card Centrifuge IIDC-2000

Blood type ID Card Centrifuge IIDC-2000 blood card centrifuge especially used for blood group testing, microcolumn gel immune detection. Spins sample volumes up to 12 and 24 cards micro gel rotor with maximum speed 3840 rpm. Driven by a brushless motor. Auto-balancing centrifuge for blood bank is particularly designed for the quick centrifuge in blood bank with storages up to 9 program records for users. For safety, we introduced an electronic door lock function and, imbalance detection, over-speed protection, inching function, short-spin function.



- Specifically for blood group test, micro-column gel immune detection
- The machine adopts brushless motor driving, micro-computer control
- Electric lid lock, super speed, and self-test protection
- Low noise, without carbon dust pollution
- Professional program design Can be run directly Without tedious parameter settings
- Maximum capacity of up to 12 and 24 cards micro gel rotor

Blood group card dedicated program automatically

Timer Range	Speed	RCF
0 to 2 min	780 rpm	75xg
2 to 3 min	1280 rpm	200xg
3 to 10 min	3840 rpm	1790xg

Blood type ID Card Centrifuge IIDC-3000

Blood type ID Card Centrifuge IIDC-3000 with standard programs are built-in for blood type determination and all kinds of experiments of serology, which standardize the test and research of clinical. Users can set the programs with a simple touch screen button. Suitable for blood group testing, microcolumn gel immune detection. Spins sample volumes up to 24 cards micro gel rotor with a maximum speed of 3000 rpm. Driven by a DC brushless motor. Auto-balancing centrifuge for blood bank is particularly designed for the quick centrifuge in blood bank with storages up to 9 program records for users. For safety, we introduced an electronic door lock function and, imbalance detection, over-speed protection, inching function, short-spin function.



- Specifically for blood group test, micro-column gel immune detection
- The machine adopts DC brushless motor driving, micro-computer control
- Electric lid lock, super speed, and self-test protection
- Low noise, without carbon dust pollution
- Professional program design Can be run directly Without tedious parameter settings
- Maximum capacity of up to 24 cards micro gel rotor

