

CENTRIFUGES



cs@irisanalytical.com



+1 706 760 4500



www.irisanalytical.com

PRP/PPP Centrifuge IPRPC-1000 designed for platelet-rich plasma (PRP) from whole blood. Used in fat transplantation surgery and autologous blood to extract for PRP cosmetic surgery. Spins sample volumes up to 6X50ml Syringe with maximum speed 4500 rpm. Driven by conversion motor. The rotors and the accessories can be disinfected under high temperature and high pressure.

- Specifically for beauty hospital
- Suitable for autologous fat transplantation surgery and autologous blood to extract for PRP cosmetic surgery
- Use of 10ml, 20ml, 50ml Syringe, and 10ml Vacutainer
- Available different PRP accessories
- Frequency conversion motor and maintenance-free
- Micro-computer control ,digital display
- Electronic lid locks and All-steel body to ensure safety
- The rotors and the accessories can be disinfected under high temperature and high pressure

