

CENTRIFUGES



cs@irisanalytical.com



+1 706 760 4500



www.irisanalytical.com



Cyto Centrifuge ICC-1000 is an independent device especially designed for cytological applications that transfer cells from fluids to glass slides for microscopic examination. The ICC-1000 has the capacity for up to 0.5 to 4 ml sample holders. Each specimen holder will accommodate a specimen of up to 12 ml. All the cells from a sample are placed on a glass slide evenly in a specific, square area; therefore, screening of the slide is shortened and easy. It is simple to load and unload chambers.

