

REFRIGERATORS & ULTRA-LOW FREEZERS



 cs@irisanalytical.com

 +1 706 760 4500

 www.irisanalytical.com



Double Door Deep Freezers -40°C IUDF-200D is a double door with dual compressors and chambers, which allows you to control the upper chamber and lower chamber autonomously. Our IUDF-200D is energy efficient, noiseless, and reliable which combines the maximum consistency and higher working performance with cost-effective operation and advanced features. The maximum capacity is up to 450L. Microcomputer temperature control system to set the inner temperature at the range of -20°C to -40°C. The built-in protective seal door to maintain a constant temperature and is easy to clean. The low temp freezer is best for the storage of various medical and scientific materials, including blood plasma, reagent, and more.

- Equipped with a broad range of capacities: 531 L
- Microprocessor-based temperature controller
- CE certified
- Temperature Range: -10 to -40°C
- Digital display with the high precision value
- Equipped with compressor, fan electromotor with EBM and unique refrigeration circulation and unipolar compressor oil-lubricated technology
- Multiple alarm and back-up battery
- Eco-friendly refrigerators: R290
- Available with the door lock, inventory boxes, insulation, shelves, ampoule, casters, compressor, fan motor
- Lower energy consumption and heat emission cost with the low noise level
- Covering with PCM and Aluminum plate with spraying and ABS drawers
- Lower energy consumption and heat emission cost with low noise level
- Standard USB for data transfer and storing information
- Build-in thermal insulating foamed door to prevent loss of refrigerating capacity in an effective way
- CFC-free polyurethane material to improve heat insulating performance to a large extent
- Compatible with the audible and visual alarm system
- Provide keyboard lock function to prevent any adjustment of operating parameters without access
- Equipped with the opaque door handle to prevent damage during transportation and save space

