

REFRIGERATORS & ULTRA-LOW FREEZERS



 cs@irisanalytical.com

 +1 706 760 4500

 www.irisanalytical.com



Ultra-Low Chest Freezers -86°C IUCF-2000 Series design for large amounts of freezer storage, the available maximum capacities are: 105 L and 485 L with a temperature range of -40°C to -85°C by offering excellent temperature uniformity and stability. Build-in inner door to avoid more cooling leaking out, and provide more stable inside temperature performance and complete temperature alarm system featuring. Used to stores vaccines, blood plasma, and many other biological products. Features with a cooling fan for instrument safety, durability, and endurance purpose. The range of medical freezers storage products will permit you to choose the perfect match to meet your requirements through cost and energy saving.

- Available in two models: single door chest and double door chest freezer
- Provide high-density CFC and Free urethane foam insulation making instrument more rigid and strong, stable storage temperature
- Temperature can be controlled using a microprocessor
- Digital display with adjustable temperature range: -40°C to -85°C
- Accommodate LED digital display for clear observation
- Build-in inner door to avoid more cooling leaking out, and provide more stable inside temperature performance
- Provide easy door opening with upright side door and handle limit door hinge to prevent a collision when door open
- Lockable caster design for customer to move as needed
- Wide voltage tolerance design with an applied voltage range of 187 to 242 V/AC
- Regulator and control password protection function to avoid unwanted access

