

REFRIGERATORS & ULTRA-LOW FREEZERS



 cs@irisanalytical.com

 +1 706 760 4500

 www.irisanalytical.com



Ultra-Low Upright Freezers -86°C IUUF-3000 suitable for blood stations, hospitals, health and epidemic prevention system, scientific research institutes, biological engineering, a laboratory in colleges and universities, military-industrial enterprises, etc. Available with maximum capacity up to 668 L. Microcomputer temperature control system to set the inner temperature at the range of -40 to -86°C.

- Equipped with a broad range of capacity: 668 L
- Microprocessor-based temperature controller
- CE certified
- Temperature Range: -40 to -86°C
- Digital display with the high precision value
- Equipped with compressor, fan electromotor with EBM and unique refrigeration circulation and unipolar compressor oil-lubricated technology
- Multiple alarm and back-up battery
- Eco-friendly refrigerators
- Available with the door lock, inventory boxes, insulation, shelves, ampoule, casters, compressor, fan motor
- Lower energy consumption and heat emission cost with the low noise level
- Covering with PCM and Aluminum plate with spraying and ABS drawers
- Lower energy consumption and heat emission cost with the low noise level
- Standard USB for data transfer and storing information
- Build-in thermal insulating foamed door to prevent loss of refrigerating capacity in an effective way
- CFC-free polyurethane material to improve heat insulating performance to a large extent
- Compatible with an audible and visual alarm system
- Provide keyboard lock function to prevent any adjustment of operating parameters without access

