

REFRIGERATORS & ULTRA-LOW FREEZERS



 cs@irisanalytical.com

 +1 706 760 4500

 www.irisanalytical.com



Blood Bank Refrigerator IUBR-2000 Series equipped with an advanced automatic defrost sensor to manage the defrost cycle with capacities of 4/66 and 12/132 blood bags with volumes of 450 ml. Intuitive two alarming of audible buzzing and visual flashing light, dual transparent insulated glass door, and display interface provide navigation to feature user interface design easily along with 5 alarm conditions including high/low temperature, power failure, voltage abnormal, door ajar, backup battery. Ergonomic design offers a clear view of stored items over large framed openings. The thin shape allows for cool recovery of your products, without the concern for swinging door clearances. Customers can choose from two types of interior shelving to suit their needs, security and safety lock to avoid unwanted access, provide storage space for easy sorting of a variety of blood products, and casters design and interior light.

- Automatic defrost maintains the optimum cooling capacity
- Designed with an acrylic inner small door,
- Door-activated ON/OFF switch and self-regulating switch on the panel also glass door with electric heat reduces condensation
- Excellent battery back up
- Positive, forced-air circulation maintains temperature uniformity throughout
- Regular cabinet temperature up to 2-6°C
- Equipped with audible buzzing and visual flashing light type alarm system
- Multiple alarm system for safe and consistent process
- Assure inner temperature uniformity by an intelligent controlled fan with forced cycle system
- Assure safety by unwanted access
- Caster design and interior light
- Excellent durability and reliability

Accessories

- 2, 4, 6 inner doors
- Testing hole
- Outdoor latch
- Light
- Thermal printer

Optional Accessories

- Remote alarm
- Recorder temperature Recorder

