

REFRIGERATORS & ULTRA-LOW FREEZERS



 cs@irisanalytical.com

 +1 706 760 4500

 www.irisanalytical.com



Double Pharmacy Refrigerator 2 to 8°C IUPR-2000D Series are 10 shelves refrigerators that are always a spectacular addition to pharmaceutical laboratories, and they are remarkably economical. This allows us a reliable and regular development of our pharmaceutical refrigerators. Each unit consists of a robust housing and quality constituents provide reliability and durability. Ergonomic design offers a clear view of stored items over large framed openings. The thin shape allows for cool recovery of your products, without the concern for swinging door clearances. Customers can choose from two types of interior shelving to suit their needs, security and safety lock to avoid unwanted access, provide storage space for easy sorting of a variety of blood products, and casters design and LED interior light.

- The cooling fan is permanently lubricated for long life
- Temperature range is up to 2 to 8°C
- Adjustable shelves and LED interior light, caster design
- Automatic defrost maintains the optimum cooling capacity
- Improved refrigeration system design for more operative cooling and quick recovery
- Designed with 1 glass door with heating reflect reduces condensation
- Excellent battery back up
- Positive, forced-air circulation maintains temperature uniformity throughout
- Temperature variation within $\pm 3^{\circ}\text{C}$
- High-efficiency field compressor with known field consistency
- Equipped with audible buzzing and visual flashing light type alarm system
- Multiple alarm system for safe and consistent process
- Assure inner temperature uniformity by an intelligent controlled fan with forced cycle system
- Assure safety by unwanted access and pharmaceutical products of many types
- Excellent durability and reliability

