

REFRIGERATORS & ULTRA-LOW FREEZERS



 cs@irisanalytical.com

 +1 706 760 4500

 www.irisanalytical.com



- Full Description
- Specification

Ultra-Low Chest Freezers -164? IUCF-5000 build-in foaming technology, and VIP insulation which greatly improves the insulation effect, VIP vacuum insulation board, tightly locks the cold air to ensure the great insulation effect. Features with single-pole lubricated compressor refrigeration technology for better-quality refrigeration capacity. Allow maximum capacity up to 128 L with a temperature range of -110°C to -164? by offering excellent temperature uniformity and stability. Available with upright racks, alarm lamp, voltage compensation, remote alarm, and liquid N2 backup system as optional accessories.

- Equipped with VIP vacuum insulation board, tightly locks the cold air
- Equipped with a broad range of capacity classes: 128 L
- Microprocessor-based temperature controller
- Temperature Range: -10 to -164°C
- Digital display
- Multiple alarm and back-up battery
- Coil evaporator
- CE certified
- Assure safety by providing keyboard lock function and password protection function
- Designed with unique secondary foaming technology, and super thick VIP insulation
- Mixture gas refrigerators
- Optional single type upright racks Optional alarm lamp, voltage compensation, remote alarm, and liquid N2 backup system

Model	IUCF-5128
Capacity	128 L
Temperature Range	-110 to -164?
Materials	Internal: structure steel, External: Steel plates with spraying
Alarms	Visual & Audio high and low-temperature alarm, High/low voltage, sensor failure alarm, Low battery, Condenser heat dissipation, System failure.
Type	Chest
Evaporator	Coil type evaporator
Display	Digital
Refrigerant	Mixture gas
Casters	6
Compressor	1 pcs
Motor	Fan Motor: EBM
Ambient Temperature	16 to 32°C
Certification	CE
Power Consumption	4600 W
Internal Dimension	490×470×585 mm
External Dimension	1510×985×1120 mm
Net Weight	305 kg
Electrical Requirement	380V/50 Hz
Catalogue No.	5001050503

