

Liquid Handling Products



cs@irisanalytical.com



+1 706 760 4500



www.irisanalytical.com



Fully Autoclavable Multichannel Micropipette AMQ Series covers all volumes from 0.5 to 300 µl. The 121°C fully autoclavable, allow easy turning and cleaning, and decrease the risk of contamination. Available with 3 different volumes ranges including 0.5 to 10 µl, 5 to 50 µl, and 50 to 300 µl. Additionally come with a finger hook take the load for extra hold, a soft pipetting button for free and comfort rotation and changing volume as well, an ejector button, digital volume display for easy readout.

- 121°C Fully autoclavable
- Dispensing head rotates for effortless pipetting convenience
- Individual piston and tip cone assembly allow for easy repair and maintenance
- Compound material-made tip cone secures high sealing performance
- Compatible with most universal tip brands
- After autoclaving, the pipette must be cooled down and left to dry for 12 hours before use.
- It is recommended to check the performance of the pipette after each autoclaving.
- Greasing and sealing of the pipette piston after every 10th autoclaving will secure an enhanced function
- Supplied with a calibration certificate of quality ISO8655
- Ergonomic, lightweight design and soft plunger movement reduce manual handling
- Fast and accurate volume setting, separate tip ejector and clear digital readout
- Provide the most accurate volume range, reliable, comfortable for extensive pipetting

Model	Volume Range(μl)	Increment(μl)	Test Volume(μl)	Accuracy		SD (μl)	CV (%)	Catalogue No.
				(?)	(%)			
AMQ-8								
AMQ-8-10	0.5-10	0.1	10	0.15	1.50	0.15	1.50	4005005-04
			5	0.125	2.50	0.125	2.50	
			1	0.04	4	0.04	4	
AMQ-8-50	5-50	0.5	50	0.5	1	0.25	0.50	4005010-04
			25	0.375	1.50	0.25	1	
			5	0.15	3	0.1	2	
AMQ-8-300	50-300	5	300	2.1	0.70	0.75	0.25	4005015-04
			150	1.5	1	0.75	0.50	
			50	0.75	1.50	0.4	0.80	
AMQ-12								
AMQ-12-10	0.5-10	0.1	10	0.15	1.50	0.15	1.50	4005020-04
			5	0.125	2.50	0.125	2.50	
			1	0.04	4	0.04	4	
AMQ-12-50	5-50	0.5	50	0.5	1	0.25	0.50	4005025-04
			25	0.375	1.50	0.25	1	
			5	0.15	3	0.1	2	
AMQ-12-300	50-300	5	300	2.1	0.70	0.75	0.25	4005030-04
			150	1.5	1	0.75	0.50	
			50	0.75	1.50	0.4	0.80	

