



Liquid Handling Products



cs@irisanalytical.com



+1 706 760 4500



www.irisanalytical.com



Electronic Micropipette EP Series is based on Stepper Motor with Microprocessor controlled piston movement. These are user friendly and offer effortless pipetting. Accuracy & Precision levels are very high and user-independent. This is a complete solution for every pipetting need.

- Pipetting modes: Standard mode, Stepper mode, Dilution mode, and Calibration mode
- Variable Working Speeds: Speed can be selected on a scale of 1-5
- Volume Selection: Easy volume selection by 'UP' & 'DOWN' keys.
- Large Fire Button: Single 'FIRE BUTTON' for filling and dispensing.
- Large Liquid Crystal Display: Two-line Alphanumeric LCD.
- Built-in streamlined tip ejector: For taken easy tip ejection and access to narrow naked bottles.
- Calibration: User-friendly calibration mode facilitates easy calibration.
- Battery: 9V Ni-MH rechargeable batteries are used for 8 hours of operation with fully charged batteries.

Model	Volume Range (μl)	Increment (μl)	Accuracy		CV		Catalogue No.
			(+%)	(±μl)	(±%)	(±μl)	
EP-10	0.2 - 10μl	0.05	1.2	0.12	0.8	0.08	4005505-04
EP-120	5 - 120μl	0.5	0.8	0.96	0.3	0.36	4005510-04
EP-300	20 - 300μl	1	0.8	2.4	0.3	0.9	4005515-04
EP-1000	50 - 1000μl	5	0.8	8	0.3	3	4005520-04
EP-5000	100 - 5000μl	10	0.8	40	0.3	15	4005525-04

