

LAMINAR FLOW CABINETS



cs@irisanalytical.com



+1 706 760 4500



www.irisanalytical.com



Horizontal Laminar Flow Cabinets LHC Series has been specially designed to provide germ-free, sterilize air conditions to the work area with lighting from a LED light included the internal brightness helps the operator.

Horizontal Laminar Flow Cabinets LHC-1000 Series

Horizontal Laminar Flow Cabinets LHC-1000 Series is a workbench, which creates a particle-free working environment by taking air through a filtration system and exhausting it across a work surface in a laminar or unidirectional air stream. The laminar flow cabinet is enclosed on the sides and kept under constant positive pressure in order to prevent the infiltration of contaminated room air. Laminar flow cabinet is widely used in medical research laboratories, hospitals, manufacturing facilities and other research and production environments. LED Display: Microprocessor control system, LED display, Waterproof socket: 2 waterproof sockets are located in the side panel, for optimum convenience of using small devices inside the cabinet, Side glass window: The transparent side glass windows maximize light and visibility inside the cabinet, providing a bright and open working environment.



Horizontal Laminar Flow Cabinets LHC-2000 Series

Horizontal Laminar Flow Cabinets LHC-2000 Series offers a particle-free working environment by taking air through a filtration system and exhausting it across a work surface in a laminar or unidirectional air stream. The laminar flow cabinet is enclosed on the sides and kept under constant positive pressure in order to prevent the infiltration of contaminated room air. LED Display: Microprocessor control system, LED display, Waterproof socket: 2 waterproof sockets are located in the side panel, for optimum convenience of using small devices inside the cabinet, UV lamp: Emission of 253.7 nanometers for most efficient decontamination, Gas Tap: Side glass window: The transparent side glass windows maximize light and visibility inside the cabinet, providing a bright and open working environment.



