

Filtration Products



cs@irisanalytical.com



+1 706 760 4500



www.irisanalytical.com



Multi-Place Vacuum Manifold IFM-1000 Series designed with a magnetic filter cup, funnel adaptors along with aluminum and long lasting stainless steel manifold, magnetic filter cup. We provide advanced Magnetic cup, one-handed vacuum filtration of liquids. The convenient capacities are 150 mL size allows for easy fit into small autoclaves and 300 ml, and 500 ml size is ideal for filtration of large samples. The magnets ensure a close-fitting, best sealing between the funnel cylinder and base, reducing the need for clamps.

- The magnetic and stainless steel funnels, the stainless steel funnels that can be flame sterilized
- Build with patented spin-lock design
- Quick, easy and safe to service, fast-assembly and disassembly
- Requires electricity for operating
- Large capacity in small foot print, exchangeable funnel adaptor design permits user to fit their desired funnel simply
- Cost-effective high capacity system
- Available capacity 150, 300 and 500 ml magnetic funnel

