

Filtration Products



cs@irisanalytical.com



+1 706 760 4500



www.irisanalytical.com

Multi-Place Vacuum Manifold IFM-2000 Series with glass filters filter cups holder assemblies for the microanalysis of liquids for biological or particulate contamination under vacuum. The 47 mm or 50 mm holder for microbiology and uses a 47 mm filter membrane. The anodized aluminum or stainless manifold secures the membrane between the funnel and glass base. The glass filter cups are of the highest borosilicate quality and resistant to temperature differences and mechanical tolerance though under high vacuum situations.

- The Borosilicate Glass and stainless steel funnels, the stainless steel funnels that can be flame sterilized
- Excellent pressure resistance and can be recycled for autoclaving under high temperature and pressure
- Build with patented spin-lock design
- Quick, easy and safe to service, fast-assembly, and disassembly
- Requires electricity for operating
- Large capacity in small footprint, exchangeable funnel adaptor design permits user to fit their desired funnel simply
- Cost-effective high capacity system
- Available capacity 150, 300 and 500 ml glass funnel

