

BATHS & CIRCULATORS



cs@irisanalytical.com



+1 706 760 4500



www.irisanalytical.com

Dry Bath Incubator DR Series

The Dry Bath Incubator is a microprocessor-controlled product using advanced thermoelectric technique. With this technique, more stable and accurate temperature control is achieved. Replacement is convenient with optional metal modules of various specifications configured for different types of test tubes.

Features for Dry Bath Incubator DR Series

- LCD display. It easy to setup and use
- Simultaneous display of set and actual time, temprature
- As it is equipped with various optional blocks
- Reliable design with inner extra temperature protection
- Heave of the block makes block change easy and convenient
- Temperature deflection adjusting meet with different demands
- Microprocessor controlled incubation time and temperature.
- Beep-signal / Stop after program completion

