

BATHS & CIRCULATORS



cs@irisanalytical.com



+1 706 760 4500



www.irisanalytical.com

The Dry Bath Incubator is a microprocessor-controlled product . With this technique, more stable and accurate temperature control is achieved.

Features

1. LED display
2. Current temperature and count down time display
3. Metal blocks free of sample counter-pollution
4. Easy replacement, cleaning and disinfecting of metal blocks
5. Built-in over-temperature protection
6. Calibrating temperature discrepancy
7. Buzzer alarm when temperature control ends
8. Only heating design and Custom blocks are available

