

PHYSICAL PROPERTIES TESTING

 cs@irisanalytical.com
 +1 706 760 4500
 www.irisanalytical.com



Digital Melting Point Apparatus MPA-200 Series ideal in the field of medicines, dry stuff, perfumes, and other organic crystal substances. The MPA-200 operated at room temperature to 320°C range a 0.1 °C display resolution as well as a variable heating rate and footing in such a way that the melting temperature can be recorded and allows you to review the process. Additionally, features including a video camera and a video replay function, FDA and CFR21 part 11 compliance. The temperature is selected, measured, and display digital making it accurate. Users can view with a magnifying lens to receive perfect observation. Designed for all laboratory needs with adequate requirements for sample processing.

- Features with PID and PWM
- LCD display
- Temperature range: RT to 400°C
- Detection of melting point of substances
- Multi-point calibration standard
- Portable and compact design
- Storage data: 1000 sets
- Support USB and RS232

