




# PHYSICAL PROPERTIES TESTING

 [cs@irisanalytical.com](mailto:cs@irisanalytical.com)  
 +1 706 760 4500  
 [www.irisanalytical.com](http://www.irisanalytical.com)



Digital Melting Point Apparatus MPA-100 Series convenient and automatic apparatus designed with safety and ease of operation in mind, making it ideal for use in education and industry alike, and touch control display on the front panel incorporates the whole operation. Features with a video camera and a video replay function which allows reviewing the melting point determination. It also suits for color samples and can able to test 4 samples simultaneously, GLP, GMP(EU, UPS & CN), and FDA 21CFR part 11 compliance.

- Designed with multi-user ID classify function
- 8-inch touch screen display
- Temperature range: 400 °C
- Manual and automatic melting point determination
- Easy to clean visible glass panel
- Storage up to 32 GB ROM
- Support USB, GLP, GMP(EU,UPS & CN) and FDA CFR21 part 11

