

LAMINAR FLOW CABINETS



cs@irisanalytical.com



+1 706 760 4500



www.irisanalytical.com



Horizontal Laminar Flow Cabinet – LHC-1000 Series is a work bench, which creates a particle-free working environment by taking air through a filtration system and exhausting it across a work surface in a laminar or unidirectional air stream. The laminar flow cabinet is enclosed on the sides and kept under constant positive pressure in order to prevent the infiltration of contaminated room air.

Laminar flow cabinet is widely used in medical research laboratories, hospitals, manufacturing facilities and other research and production environments.

LED Display: Microprocessor control system, LED display

Waterproof socket: 2 waterproof sockets are located in the side panel, for optimum convenience of using small devices inside the cabinet

Side glass window: The transparent side glass windows maximize light and visibility inside the cabinet, providing a bright and open working environment.

