

LAMINAR FLOW CABINETS



cs@irisanalytical.com



+1 706 760 4500



www.irisanalytical.com



Horizontal Laminar Flow Cabinet LHS-1000 Series offers particle-free working environment by taking air through a filtration system and exhausting it across a work surface in a laminar or unidirectional air stream. The laminar flow cabinet is enclosed on the sides and kept under constant positive pressure in order to prevent the infiltration of contaminated room air.

LED Display: Microprocessor control system, LED display

Waterproof socket: 2 waterproof sockets are located in the side panel, for optimum convenience of using small devices inside the cabinet

UV lamp: Emission of 253.7 nanometers for most efficient decontamination

Gas Tap:

Side glass window: The transparent side glass windows maximize light and visibility inside the cabinet, providing a bright and open working environment.

