



DRY BLOCK HEATER/THERMO MIXERS



cs@irisanalytical.com



+1 706 760 4500



www.irisanalytical.com

DR-2000 Programmable dry bath incubator applies the thermal cycle's heating lid technology, Unique dual block design operates via touch screen which is the most advanced and luxury dry block heating product in the world.

Features:

- Touch screen allows for easy to read and easy programming and easy to use
- Hold a heating lid, 60 user programmable settings
- Programmable system, program link functions, it can be runned some simple PCR tests
- Block extraction tool is supplied, allowing blocks to be removed easily

Model	DR-2000	
Temperature Control Range	RT+5°C to 100°C	
Timing Range	1min to 99h59min	
Uniformity: within one block @ 37°C	±0.2°C	
Across similar blocks@ 37°C	±0.3°C	
Temp. Control Accuracy	±0.5°C	
Display Accuracy	0.1°C	
Heating Time (25°C to Max. Temp.)	?15min	
Dimension	420x225x145 mm	
Net Weight	5.5kg	
Electrical Requirement	11V/60Hz	220V/50Hz
Catalogue No.	1603505-01	1603505-02

