

Liquid Handling Products



cs@irisanalytical.com



+1 706 760 4500



www.irisanalytical.com



Fully Autoclavable Fixed Volume Micropipette AFG Series are designed for many different types of routine laboratory work. Build-in ergonomic portable fits comfortably in the user's hand. The basic features of fixed volume pipettes are ideal for clinical diagnostics, control analysis, etc. The AFG Series Fixed Volume "Fully Autoclavable" Micropipettes are available in all capacities from 0.5 µl to 10000 µl.

Soft Grip & Ergonomic Design: New body design with soft grip and improved ergonomic for more comfort and less fatigue during operation.

Large Grippy: Specially designed large grippy is provided for a comfortable grip and ease during working.

Tip Ejector: Built-in streamlined tip ejector facilitates easy tip ejection and access to narrow-necked bottles and tubes.

Pipette Holder: The holder provided with the pipette allows easy, efficient, and safe storage.

Tip Ejection Knob: Placement and good design allow easy single-handed tip ejection.

Colour Coding: Attractive color coding for easy identification between different volume ranges.

Model	Volume (µl)	Accuracy		CV		Catalogue No.
		±%	±µl	±%	±µl	
AFG-1	1	5	0.1	5	0.1	4001005-04
AFG-2	2	4	0.1	2	0	4001010-04
AFG-5	5	2	0.1	1	0.1	4001015-04
AFG-10	10	1	0.1	1	0.1	4001020-04
AFG-20	20	1	0.2	0	0.1	4001025-04
AFG-25	25	1	0.2	0	0.1	4001030-04
AFG-50	50	1	0.4	0	0.2	4001035-04
AFG-100	100	1	0.6	0	0.2	4001040-04
AFG-200	200	1	1.2	0	0.4	4001045-04
AFG-250	250	1	1.5	0	0.5	4001050-04
AFG-500	500	1	3	0	1	4001055-04
AFG-1000	1000	1	6	0	2	4001060-04
AFG-2000	2000	1	12	0	4	4001065-04
AFG-5000	5000	1	30	0	10	4001070-04
AFG-10000	10000	0	40	0	20	4001075-04

AF Series Fixed Volume Micropipette - Fully Autoclavable

AF Series Fixed Volume "Fully Autoclavable" Micropipettes are designed for many different types of routine laboratory work. Using the same state-of-art design and basic features of Variable Volume Pipettes, these fixed volume pipettes are ideal for Clinical Diagnostics, Control Analysis etc.



The AF Series Fixed Volume "Fully Autoclavable" Micropipettes are available in all capacities from 0.5 μ l to 10000 μ l (10 ml).

Soft Grip & Ergonomic Design: New body design with soft grip and improved ergonomic for more comfort and less fatigue during operation.

Large Grippy: Specially designed large grippy is provided for comfortable grip and ease during working.

Tip Ejector: Built-in streamlined tip ejector facilitates easy tip ejection and access to narrow necked bottles and tubes.

Pipette Holder: The holder provided with the pipette allows easy, efficient and safe storage.

Tip Ejection Knob: Placement and good design allows easy single handed tip ejection.

Colour Coding: Attractive colour coding for easy identification between different volume ranges.

