

Liquid Handling Products



cs@irisanalytical.com



+1 706 760 4500



www.irisanalytical.com



Fully Autoclavable Multichannel Micropipette AMP Series has a unique structure, certifying easy tip mounting and discharging. Both the 8- and 12-channel pipettes have an external discharging knob; deliver ergonomic and easy tip ejection. This unique construction of AMP multichannel pipettes enables the thumb to rest under tip ejection, significantly decrease the thumb load. The AMP series Multichannel Pipettes are designed especially for ELISA, Molecular screening, Kinetic studies, DNA amplifications etc

- Easy Volume Adjustment: Easy volume adjustment by simply turning the plunger. the plunger does not snag gloves.
- Volume Setting: Soft click sound at every increment ensures perfect volume setting and no accidental volume change. It also facilitates single-handed operation.
- Soft Grip & Ergonomic Design: New body design with soft grip and improved ergonomics for more comfort and less fatigue during operation.
- Digital Display: 4 digit display with subdivisions provides small volume increments.
- Large Grippy: Specially designed large grippy is provided for a comfortable grip and ease during working.
- Tip Ejector: Built-in streamlined tip ejector facilitates easy tip ejection and access to narrow-necked bottles and tubes.
- User Adjustment: Re-calibration can be performed easily by the user
- Pipette Holder: The holder provided with the pipette allows easy, efficient, and safe storage.
- Easy Tip Ejection: Specially designed sequential tip ejector allows effortless tip ejection.
- Operation Flexibility: Lower housing can be rotated 360° providing flexibility during operation.
- Tip Ejection Knob: Placement and good design allow easy single-handed tip ejection.
- Colour Coding: Attractive color coding for easy identification between different volume ranges.

8 - Channel Fully Autoclavable Micropipette							
Model	Volume Range (µl)	Increment(µl)	Accuracy		CV		Catalogue No.
			±%	±µl	±%	±µl	
AMP-8-10	0.5 - 10µl	0.02	1.6	0.16	1	0.1	4004505-04
AMP-8-20	02 - 20µl	0.02	0.8	0.16	0.4	0.08	4004510-04
AMP-8-50	05 - 50µl	0.1	0.8	0.4	0.4	0.2	4004515-04
AMP-8-100	10 - 100µl	0.2	0.8	0.8	0.3	0.3	4004520-04
AMP-8-200	20 - 200µl	0.2	0.8	1.6	0.3	0.6	4004525-04
AMP-8-300	40 - 300µl	0.2	0.8	2.4	0.3	0.9	4004530-04
12 - Channel Fully Autoclavable Micropipette							
AMP-12-10	0.5 - 10µl	0.02	1.6	0.16	1	0.1	4004535-04
AMP-12-20	02 - 20µl	0.02	0.8	0.16	0.4	0.08	4004540-04
AMP-12-50	05 - 50µl	0.1	0.8	0.4	0.4	0.2	4004545-04
AMP-12-100	10 - 100µl	0.1	0.8	0.8	0.3	0.3	4004550-04
AMP-12-200	20 - 200µl	0.1	0.8	1.6	0.3	0.6	4004555-04
AMP-12-300	40 - 300µl	0.1	0.8	2.4	0.3	0.9	4004560-04

