



Liquid Handling Products



cs@irisanalytical.com



+1 706 760 4500



www.irisanalytical.com



- Motor Controlled piston movement.
- Re-circulation system allows quick purging without loss of reagent.
- Easy, user-friendly, TFT touchscreen guides you to perform quick and accurate titration.
- Three dispensing speeds including drop-wise dispensing allows you to achieve the end point very accurately.
- Dispensing can be performed by a touch of the screen and stops immediately when the touch is removed.
- 20 results can be saved in memory.
- Computer connectivity.
- Conforms to ISO 8655 standards.

Component	Description
Piston	PTFE & ETP
Cylinder	Borosilicate Glass
Valve Housing	PFA
Re-Circulation Valve Housing	PFA
Valve Assembly	Borosilicate Glass Ball & Seat
Discharge Assembly	PTFE
Delivery Tube	FEP
Inlet Tube	FEP
Calibration	Individually calibrated and certified.
Accuracy & Reproducibility	In accordance with ISO 8655 standards

