

Liquid Handling Products



cs@irisanalytical.com



+1 706 760 4500



www.irisanalytical.com



SFP Series Single Channel Fixed Volume Pipettors features twelve fixed volume micropipettes. These pipettes are calibrated in accordance ISO/DIS 8655. Each pipettor passes a quality control/calibration test and is supplied with a certificate.

- Ergonomic design provides excellent operating experience
- Large display window allows for easy volume identification
- Autoclavable tip-cone
- Integral tip ejector for access into narrow tubes
- Easy to calibrate and maintain with tool supplied
- Used with 'universal' pipette tips
- Supplied with an individual calibration certificate

All Iris pipettes have been quality tested according to ISO8655-2:2002 and are supplied with individual calibration certificates. The quality control includes gravimetric testing of each pipette with distilled water at 22°C.

