



Liquid Handling Products



cs@irisanalytical.com



+1 706 760 4500



www.irisanalytical.com





AVP Series Variable Volume Fully Autoclavable Micropipettes are available in eight variable volume from 0.1µl to 10ml. The volume adjustment is selected by turning the thumb plunger mechanism. Iris pipettes have been quality tested according to ISO8655-2:2002 and are supplied with individual calibration certificates. The quality control includes gravimetric testing of each pipette with distilled water at 22°C.

The Iris AVP pipettes are fully autoclavable, enable easy cleaning and reduce the risk of contamination. Steam autoclaving can be performed at 121°C, 1 bar for 20 minutes. After autoclaving, the pipette must be cooled down and left to dry for 12 hours before use.

- Fully autoclavable at 121°C
- Pipettes cover a volume range from 0.1µL to 10mL
- Ergonomic design provides excellent operating experience
- Large display window allows for easy volume identification
- Integral tip ejector for access into narrow tubes
- Easy to calibrate and maintain with tool supplied
- Autoclavable tip-cone
- Used with 'universal' pipette tips
- Each pipettes supplied with an individual calibration certificate

