



PUMPS & ASPIRATORS



cs@irisanalytical.com



+1 706 760 4500



www.irisanalytical.com



Compact Peristaltic Pumps are equipped with an LCD display controller that permits monitoring of parameters and changes of speed. Appropriate for continuous dosing application. A quick exchange of pump heads without the use of tools.

Iris Compact Peristaltic Pump ICPP-1000 Series is perfect for continuous dosing applications. Unique LED touch screen display function as well as digital knob and mechanical control mode and 3 pump head drive of 0.1 to 350 rpm speed range.

Compact Peristaltic Pump ICPP-1000 Series

Compact Peristaltic Pump ICPP-1000 Series is perfect for continuous dosing applications. Unique OLED touch screen display function as well as digital knob and mechanical control mode and 3 pump head drive of 0.1 to 350 rpm speed range.



Compact Peristaltic Pump ICPP-1445-1PH has 0.004 to 149.23 ml/min flow rates and excellent dosing that used for timer function of simple dispensing function along with external control 0 to 5V, 4 to 20 mA, and 0 to 10V (OPTIONAL) USB connector, RS485 interface.

Compact Peristaltic Pump ICPP-1445-2PH has 0.0024 to 930 ml/min flow rates and excellent dosing that used for timer function of simple dispensing function along with external control 0 to 5V, 4 to 20 mA, and 0 to 10V (OPTIONAL) USB connector, RS485 interface.

Compact Peristaltic Pump ICPP-1445-1PH

Compact Peristaltic Pump ICPP-1445-1PH has 0.004 to 149.23 ml/min flow rates and excellent dosing that used for timer function of simple dispensing function along with external control 0 to 5V, 4 to 20 mA, and 0 to 10V (OPTIONAL) USB connector, RS485 interface.



- Low power consumption
- LED display the current speed
- ABS plastics constructing material
- A quick exchange of pump heads without the use of tools

Compact Peristaltic Pump ICPP-1445-2PH

Compact Peristaltic Pump ICPP-1445-2PH has 0.0024 to 930 ml/min flow rates and excellent dosing that used for timer function of simple dispensing function along with external control 0 to 5V, 4 to 20 mA, and 0 to 10V (OPTIONAL) USB connector, RS485 interface.



- Low power consumption
- LED display the current speed
- ABS plastics constructing material
- A quick exchange of pump heads without the use of tools

