



PUMPS & ASPIRATORS



cs@irisanalytical.com



+1 706 760 4500



www.irisanalytical.com

Standard Peristaltic Pump ISPP-3888P has 0.000829 to 2280 ml/min flow rates and excellent dosing that used for the timer function of simple dispensing function along with 8 modes of external control of RS232 standard interface along with optional RS485 Modbus interface communicator.

- Support RS232 and RS485, the pump can be connected with PC to reach dispensing and filling function
- LED display the current speed
- Plastic coated metal housing constructing material
- Mechanical keypad to control the start/stop, speed, and direction
- Negligible noise and vibration
- Design with ultra-quiet drive and precision control
- Testing time range 0.5 s to 999 s
- Left-right reversal
- Display lamp of external mode
- Equipped with timer for 60 seconds test
- A quick exchange of pump heads without the use of tools

