



PUMPS & ASPIRATORS



cs@irisanalytical.com



+1 706 760 4500



www.irisanalytical.com

Medium Flow Rate Industrial Peristaltic Pump IFPP-220 has 0.000166 to 1330 ml/min flow rates and excellent dosing that used for timer function of simple dispensing function along with RS232, RS485 Modbus interface communicator.

- Allow stacking of 4 pumps head for same drive, creating multichannel pump system
- 4.3-inch industrial grade LCD screen display
- High repeatability on all channels
- The directional signal comes with passive switch signal like foot pedal along with optional active switch signal: 5V, 12V, 24V
- Flow rate accuracy of pump $<\pm 0.5\%$ and 300,000 times keypad lifetimes
- To control the speed of the pump knob are used
- Equipped with 3 different working mode including fixed volume metering, fixed time and volume, and timer start and stop
- Terrific online micro-adjusting and intelligent calibration function
- Easy installation and low maintenance condition
- Materials have good rigidity and high precision along with resisting organic solvent corrosion property
- Adjust the pump as per the need

