



PUMPS & ASPIRATORS



cs@irisanalytical.com



+1 706 760 4500



www.irisanalytical.com

Low Peristaltic Pump IDPP-590 has 0.00166 to 2280 ml/min flow rates and excellent dosing that used for the timer function of simple dispensing function along with RS232, RS485 Modbus interface communicator.

- Used to pre-set dispensing volume, dispensing time, pause time, copy numbers and back suction angle
- Fabricated with a stepper motor that runs specific control
- Speed range is 1-600 rpm
- An automatic temperature control system
- Design with 4 pump head with excellent flow rate and accurate dispensing volume range
- Excellent anti-corrosive and anti-static function

