

Tabletop Autoclave ITA-314



The ITA-314 front loading 'Benchtop' autoclave range is available in 14 L size. These units require no plumbing, and only a standard single phase power supply. Build in LCD display with 8 programs automatically feeding water and draining water system which is easy to operate and complete automatic operation, check and monitor by computer.

Technical Specification:

Model	ITA-314
Capacity (liters)	14 L
Chamber Dimensions (mm)	250×270 mm
Vacuum Drying	Yes
External Dimensions (w x d x h) mm	480×480×420 mm
Net Weight (kg)	42.5 kg
Gross Weight (kg)	50.5 kg
Shipping Dimensions (w x d x h) mm	630×600×500 mm
Power (W)	1600 W
Catalogue No.	5801005-04

Technical Data:

Model	ITA-314	ITA-317	ITA -323
Chamber Capacity	14 L	17 L	23 L
Sterilizer Construction	121/134		
Stainless steel instrument container	Stainless steel SUS-304		
Dry Procedure	Vacuum drying		
Minimum Desk Dimension	480×350 mm	480×400 mm	480×500 mm
Heater	1600 W	1800 W	2200 W
Net Weight	42.5 kg	47 kg	53 kg
Gross Weight	50.5 kg	55 kg	62 kg
Chamber Dimensions	250×270 mm	250×350 mm	250×490 mm
Packing Dimensions	630×600×500 mm	630×600×500 mm	830×600×500 mm
Overall Dimensions	480×480×420 mm	480×530×420 mm	480×690×420 mm
Electrical Requirements	-	-	-
Catalogue No.	5801005-04	5801010-04	5801015-04



Iris Analytical Ltd.

66 W Flagler St #900 Miami, Florida 33130 United States