

Laboratory Glassware Washers

IGW-2200(380V/50Hz)



Laboratory Glassware Washers IGW-2200 is the most popular glassware washer that provides washing and drying function with HEPA filter can clean and dry the inner of tubular things especially for laboratory and medical applications.

The IGW-2200 is decontaminating equipment that worked on moist heat to remove contamination and bacteria with a large wash cabinet and short programme cycles. The key features include stainless steel interior, exterior, racks and inserts withstand heavy use and provide protection against corrosion resistance materials and offer internal components peristaltic pumps and seals laboratory-grade for excellent durability of instruments.

Technical Specification:

Model	IGW-2200
Functions	Washing
Capacity	198 L
Materials	Internal Panel: 316L, External Panel: 304L
Peristaltic pump	2 to 4 pcs
Programs	36 Default, 100 Custom
Interface	RS232, pure booster pump
Heating Power	3/9/18 KW
Dimensions	1000×617×770 mm
Electrical Requirements	380V/50Hz
Catalogue No.	370100503

Technical Data:

Model	IGW-2200	IGW-2400
Function	Washing	Washing + Drying with HEPA filter
Capacity	198 L	198 L
Washing Temperature	? 93°C	? 95°C
Materials	Internal Panel: 316L, External Panel: 304L	Internal Panel: 316L, External Panel: 304L
Peristaltic pump	2 to 4 pcs	2 to 4 pcs
Sensor	PT1000	PT1000
Cleaning Quantity	Sampling Tube: 448 or 476 pcs, Measuring Tube 100 (ml): 84 pcs, Pipette: 238 pcs	Sampling Tube: 448 or 476 pcs, Measuring Tube 100 (ml): 84 pcs, Pipette: 238 pcs
Programs	36 Default, 100 Custom	36 Default, 100+Custom
Interface	RS232, pure booster pump	RS232, pure booster pump
Heating Power	3/9/18 KW	3/8/18 KW
Total Power	5/11/20 KW	5/11/20 KW
Dimension	1000×617×770 mm	1000×935×770 mm
Electrical Requirement	220V/50Hz	220V/50Hz
Catalogue No.	370100502	370101002



Iris Analytical Ltd.

66 W Flagler St #900 Miami, Florida 33130 United States