

Laboratory Glassware Washers

IGW-3400(220V/50Hz)



Laboratory Glassware Washers IGW-3400 designed to provide excellent cleaning and drying efficiency with HEPA filter, can clean and dry the inner of tubular things especially for laboratory and medical applications.

The IGW-3400 is decontaminating equipment that worked on moist heat to remove contamination and bacteria with a large wash cabinet and short programme cycles. The key features include stainless steel interior, exterior, racks and inserts withstand heavy use and provide protection against corrosion resistance materials and offer internal components peristaltic pumps and seals laboratory-grade for excellent durability of instruments.

Technical Specification:

Model	IGW-3400
Functions	Washing + Drying with HEPA filter
Capacity	160 L
Circulating System	0 to 400 L/min
Materials	Internal Panel: 316L, External Panel: 304L
Peristaltic pump	2 to 3 pcs
Programs	36 Default, 50+user defined
Interface	RS232
Heating Power	3 KW/9 KW
Dimensions	835×617×765 mm
Electrical Requirements	220V/50Hz
Catalogue No.	370151002

Technical Data:

Model	IGW-3200		IGW-3400	
Function	Washing		Washing + Drying with HEPA filter	
Capacity	126 L		160 L	
Circulating System	0 to 400 L/min		0 to 400 L/min	
Washing Temperature	? 93°C		? 93°C	
Materials	Internal Panel: 316L, External Panel: 304L		Internal Panel: 316L, External Panel: 304L	
Peristaltic pump	2 pcs		2 to 3 pcs	
Sensor	PT1000			
Programs	36 Default, 50+user defined		36 Default, 50+user defined	
Interface			RS232	
Installation	Counter top		Under counter	
Heating Power	2 KW		3 KW/9 KW	
Total Power	3 KW		4/10 KW	
Dimension	680×617×765 mm		835×617×765 mm	
Electrical Requirement	220V/50Hz	380V/50Hz	220V/50Hz	380V/50Hz
Catalogue No.	370150502	370150503	370151002	370151003



Iris Analytical Ltd.

66 W Flagler St #900 Miami, Florida 33130 United States