

Laboratory Glassware Washers

IGW-4000(220V/50Hz)



Laboratory Glassware Washers IGW-4000 closed and compact sophisticated model and easy to use. The chamber capacity is up to 308 L and perfect solution for both small and large items as it has countertop installation as well as under counter installation.

The IGW-4000 is decontaminating equipment that worked on moist heat to remove contamination and bacteria with a large wash cabinet and short programme cycles. The key features include stainless steel interior, exterior, racks and inserts withstand heavy use and provide protection against corrosion resistance materials and offer internal components peristaltic pumps and seals laboratory-grade for excellent durability of instruments.

Technical Specification:

Model	IGW-4000
Capacity	308 L
Materials	Internal Panel: 316L, External Panel: 304L
Programs	36 Default, 100+Custom
Interface	RS232, pure booster pump
Heating Power	3/8/18 KW
Dimensions	1350×925×760 mm
Electrical Requirements	220V/50Hz
Catalogue No.	370200502

Technical Data:

Model	IGW-4000	
Capacity	308 L	
Materials	Internal Panel: 316L, External Panel: 304L	
Cleaning Quantity	Sampling Tube: 672 pcs, Measuring Tube 100 (ml): 126 pcs, Pipette: 238 pcs	
Programs	36 Default, 100+Custom	
Interface	RS232, pure booster pump	
Heating Power	3/8/18 KW	
Total Power	5/11/20 KW	
Dimension	1350×925×760 mm	
Weight		
Electrical Requirement	220V/50Hz	380V/50Hz
Catalogue No.	370200502	370200503



Iris Analytical Ltd.

66 W Flagler St #900 Miami, Florida 33130 United States