

# Portable Centrifuge IPC-130(220V/50Hz)



Portable Centrifuge IPC-130 s are ideal for the separation of serum plasma, urea, blood samples and other routine applications in hospital and research laboratories.

The IPC-130 are most powerful and competitive centrifuge in its class. it can be placed with 6 of 15ml conical tubes or glass tubes, 8 of 15ml conical tubes or glass tubes conical tubes which with screw caps The speed can be 100- 5000RPM. the LCD display with ability to change be between RPM and RCF.

## Technical Specification:

|                         |                     |
|-------------------------|---------------------|
| Product Type            | Portable Centrifuge |
| Model                   | IPC-130             |
| Speed Range (rpm)       | 5000 rpm            |
| Max RCF                 | 3074xg              |
| Capacity                | 5mlX12              |
| Power                   | 40W                 |
| Electrical Requirements | 220V/50Hz           |
| Catalogue No.           | 140201502           |

## Technical Data:

| Model                                    | IPC-110       |           | IPC-120   |           | IPC-130   |           |
|--|---------------|-----------|-----------|-----------|-----------|-----------|
| Max. capacity                            | 15mlX6        |           | 15mlX8    |           | 5mlX12    |           |
| Max. speed rpm (Angle rotor)             | 5000 rpm      |           |           |           |           |           |
| Speed increment                          | 100 rpm       |           |           |           |           |           |
| Max. RCF[xg]                             | 3074xg        |           |           |           |           |           |
| Noise level at max. speed(approximately) | ?60dB(A)      |           |           |           |           |           |
| Max. Acceleration time                   | 20-80 seconds |           |           |           |           |           |
| Max. Deceleration time                   | 20-80 seconds |           |           |           |           |           |
| Power consumption drive                  | 40W           |           |           |           |           |           |
| Electrical Requirements                  | 110V/60Hz     | 220V/50Hz | 110V/60Hz | 220V/50Hz | 110V/60Hz | 220V/50Hz |
| Catalogue No.                            | 140200501     | 140200502 | 140201001 | 140201002 | 140201501 | 140201502 |



**Iris Analytical Ltd.**

66 W Flagler St #900 Miami, Florida 33130 United States