

Portable Centrifuge IPC-500



Portable Centrifuge IPC-500 low speed, user friendly large digital display to display real time information and parameter. Spins samples volume up to 12x10ml with maximum speed 4000 rpm and available with acceleration and for deceleration profile with high efficiency. Replacing and changing rotor by simply lift the rotor out of the shaft without breaking your back or damaging the centrifuge. Ideal for microfiltration and rapid spinning accelerates and decelerates within couples of minutes to ensure excellent consequences of sample separation.

Technical Specification:

Product Type	Portable Centrifuge
Model	IPC-500
Speed Range (rpm)	4000 rpm
Max RCF	1880xg
Capacity	12x10ml
External Dimensons (HxLxW) mm	300x260x210 mm
Net Weight (kg)	7 kg
Gross Weight (kg)	10 kg
Shipping Dimensions (HxLxW) mm	380x350x300 mm
Power	80 W
Electrical Requirements	AC 220V±22V□50/60Hz□5A
Catalogue No.	140400502

Technical Data:

Model	IPC-500
Capacity	12x10ml
Max RCF	1880xg
Speed Range	4000rpm
Noise	? 65dB?A?
Rotor Diameter	280 mm
Timer	1 to 99min
Power Consumption	80 W
Dimension	300x260x210 mm
Packing Dimension	380x350x300 mm
Net Weight	7 kg
Gross Weight	10 kg
Electrical Requirements	AC 220V±22V?50/60Hz?5A
Catalogue No.	140400502

Accessories:

Model	Image	Name	Details
1404005-04		Fixed Angle 1	8x15ml; 4000rpm/1880xg
1404005-04		Fixed Angle 2	12x10ml; 4000rpm/1880xg



Iris Analytical Ltd.

66 W Flagler St #900 Miami, Florida 33130 United States