

Low Speed Centrifuges IBLC-2300



Low Speed Centrifuges IBLC-2300 accommodates angle, swing rotor and swing-out ELISA plate rotor as these motor capable of accommodating a great number of tubes at once. Suitable for clinical medicine, radio-immunity, biochemistry, food hygiene, environmental protection and other research areas. Available with 9 profiles for acceleration and for deceleration with high efficiency. Features with digital LCD screen to show various parameters including running time, speed, RCF.

Rotors

Rotor	Capacity	Max Speed	Max RCF
NO. 1 Swing-out Rotor	4x100ml, 4x50ml	5000 rpm	4650Xg
NO. 2 Swing-out Rotor	8x50ml, 8x100ml, 24x15ml, 32x15ml	4000 rpm	
NO. 3 Swing-out Rotor	4x250ml (with a wind cover)		4000 rpm 3100Xg
	Adaptor	4x50 ml, 8x50 ml(Round bottom), 4x100 ml,36x10ml, 40x7 ml, Blood collection tube	
	Special hang rack	48x7 ml Blood collection tube, 64x7 ml Blood collection tube, 8x50 ml	
NO. 4 Swing-out ELISA plate Rotor	2x2x96 holes	4000rpm	1970Xg
NO.5 Swing-out Rotor	Round hang cup	4x500ml	4000 rpm 3400Xg
	Adaptor	4x250ml, 12x50ml, 36x15ml, 76x7ml Blood collection tube	
	Square hang basket	16x50ml, 40x15ml,80x7ml Blood collection tube, 112x7ml Blood collection tube, 100x1.5ml	
	Hang cup	148x5ml Radio-immunity tube, 96x7ml Blood collection tube	
	ELISA plate rack	4x2x96 holes	
NO. 6 Angle Rotor(IBLC-1300)	12x15ml	6000 rpm	5120Xg

Technical Specification:

Model	IBLC-2300
Speed Range (rpm)	6000 rpm
Max RCF	5200Xg
Capacity	4x500ml(4000rpm)
Motor Type	frequency conversion motor
Timer	1 to 99 min
External Dimemsions (HxLxW) mm	550x450x365 mm
Net Weight (kg)	53 kg
Electrical Requirements	AC 220V 50HZ 10A
Catalogue No.	140951502

Technical Data:

Model	IBLC-2100	IBLC-2200	IBLC-2300
Capacity	4X250ml(4000rpm)	4x500ml(4000rpm)	
Max RCF	4650Xg	5200Xg	
Speed Range	5000 rpm	5000 rpm	6000 rpm
Noise	?60dB?A?		
Revolutions/min	±20 rpm		
Motor	frequency conversion motor		
Timer	1 to 99 min		
Dimension	550x450x365 mm	550x450x365 mm	
Weight	40 kg	53 kg	
Electrical Requirement	AC 220V 50HZ 10A		
Catalogue No.	140950502	140951002	140951502



Iris Analytical Ltd.

66 W Flagler St #900 Miami, Florida 33130 United States